

ПЕРСОНАЛЬНАЯ ЗАЩИТА



RU VM® Перчатки смотровые нитриловые

- Представляют собой нестерильные медицинские перчатки для проведения медицинских осмотров, диагностических и терапевтических процедур, а также для обращения с зараженными медицинскими материалами

Материалы:
БНК. Без латекса

Стерилизация:
нестерильные

Тип/длина (мин. мм):
Стандарт: 240
Премиум: 240
Премиум плюс: 290
Особопрочные удлиненные (повышенного риска): 290

Внешняя поверхность:
текстурированная

Внутренняя поверхность:
хлорированная;
неопудренная

Манжета:
окантовка

Форма:
на обе руки

Количество штук в боксе/коробе:
100 / 1000; XL: 90/900 (кроме High-Risk)

Соответствие стандартам:
EN 455, ASTM D6319, ISO 11193-1

Класс риска:
Is

Упаковка:
бокс / картон

Срок годности:
5 лет

EN VM® Nitrile Examination Gloves

- Are non-sterile medical gloves, intended for conducting medical checkups, diagnostic and therapeutic procedures and for handling contaminated medical material
- Ideal for medical workers allergic to latex
- The „High-Risk“ examination gloves are designed for manipulations in condition of increased risk of infection

Materials:
NBR
Latexfree

Sterilization:
Non-sterile

Type/length (min. mm):
Standard: 240
Premium: 240
Premium Plus: 290
High-Risk: 290

External surface:
Textured

Internal surface:
Chlorinated, powder free

Cuff:
Beaded cuff

Shape:
Ambidextrous fitting

Quantity of pieces in box/box:
100 / 1000; XL: 90/900 (except High-Risk)

Compliance standards:
EN 455, ASTM D6319,
ISO 11193-1

Risk classification:
Is

Packaging:
Box / carton

Shelf life:
5 years

VM® Перчатки смотровые нитриловые

Каталожный номер	Размер	Размеры перчаток (мин. мм)		Толщина стенки (мин. мм)		Цвет	Кол-во бокс / короб
		Длина	Ширина	Палец	Ладонь		
Стандарт ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные							
1314320	XS	230	≤ 80	0,05	0,05	голубой	100/1000
1314321	S	230	80 ± 10	0,05	0,05	голубой	100/1000
1314322	M	230	95 ± 10	0,05	0,05	голубой	100/1000
1314323	L	230	110 ± 10	0,05	0,05	голубой	100/1000
1314324	XL	230	≥ 110	0,05	0,05	голубой	90/900
Стандарт ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные							
1314320V	XS	240	≤ 80	0,05	0,05	фиолетовый	100/1000
1314321V	S	240	80 ± 10	0,05	0,05	фиолетовый	100/1000
1314322V	M	240	95 ± 10	0,05	0,05	фиолетовый	100/1000
1314323V	L	240	110 ± 10	0,05	0,05	фиолетовый	100/1000
1314324V	XL	240	≥ 110	0,05	0,05	фиолетовый	90/900
Стандарт ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные № 200							
1314320Q200	XS	230	≤ 80	0,05	0,05	темно-голубой	200/1000
1314321Q200	S	230	80 ± 10	0,05	0,05	темно-голубой	200/1000
1314322Q200	M	230	95 ± 10	0,05	0,05	темно-голубой	200/1000
1314323Q200	L	230	110 ± 10	0,05	0,05	темно-голубой	200/1000
1314324Q200	XL	230	≥ 110	0,05	0,05	темно-голубой	180/1800
Стандарт ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные ♦ ЭКОФРОЙНД (биоразлагаемые)							
1314320D	XS	240	≤ 80	0,07	0,05	пастельно-зеленый	100/1000
1314321D	S	240	80 ± 10	0,07	0,05	пастельно-зеленый	100/1000
1314322D	M	240	95 ± 10	0,07	0,05	пастельно-зеленый	100/1000
1314323D	L	240	110 ± 10	0,07	0,05	пастельно-зеленый	100/1000
1314324D	XL	240	≥ 110	0,07	0,05	пастельно-зеленый	90/900
Премиум ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные							
1314325	XS	240	≤ 80	0,11	0,08	голубой	100/1000
1314326	S	240	80 ± 10	0,11	0,08	голубой	100/1000
1314327	M	240	95 ± 10	0,11	0,08	голубой	100/1000
1314328	L	240	110 ± 10	0,11	0,08	голубой	100/1000
1314329	XL	240	≥ 110	0,11	0,08	голубой	90/900
Премиум Плюс ♦ Неопудренные ♦ Текстурированные ♦ Нестерильные							
1314330	XS	280	≤ 80	0,12	0,09	голубой	100/1000
1314331	S	280	80 ± 10	0,12	0,09	голубой	100/1000
1314332	M	280	95 ± 10	0,12	0,09	голубой	100/1000
1314333	L	280	110 ± 10	0,12	0,09	голубой	100/1000
1314334	XL	280	≥ 110	0,12	0,09	голубой	90/900