



VM® Injection needles

The Vogt Medical injection needles are produced according to ISO standards and are applicable for all indications. Needles are available in all standard sizes and are marked according to the internationally accepted colour code. The triple bevel of the needle enables smooth tissue penetration and reduces tissue damage.

Description:

- Injection needles are used together with syringes, transfusion and infusion sets for the injection of medications, for taking blood and for blood transfusion
- The triple bevel and the highly polished surface of the needle enables smooth tissue penetration and reduces tissue damage
- The range of needle tip bevels (regular, short, intradermal) enable selection of the needle for every treatment according to the needs of the patient:
 - Regular bevel for intramuscular and subcutaneous injections
 - Short bevel for intravenous and intra-arterial injections and for peripheral nerve block
 - Intradermal bevel for injection of medications between the epidermis and the dermis
- Colour-coded hub for easy identification of needle size
- For use with Luer Slip and Luer Lock hubs

Medical grade materials (safety tested):

- Hub and protective cap – polypropylene
- Needle – high-quality steel with high chrome-nickel content, manufactured in compliance with AISI standards
- Needle/hub bonding – epoxy resin

Packaging:

- Blister/box/carton

Quality:

- Conforms to ISO 7864, ISO 6009 standards

Sterilisation:

- Ethylene oxide

Shelf life:

- 5 years





VM® Injection needles

Catalog No.	Gauge	Length (inch)	Outer diameter (mm)	Length (mm)	Wall	Bevel	Colour code	Quantity box/carton
1310320	20G	1 1/2"	0.9	40	thin	regular	yellow	100/5000
1310321	21G	1 1/2"	0.8	40	thin	regular	dark-green	100/5000
1310322	22G	1 1/4"	0.7	30	thin	regular	black	100/5000
1310323	23G	1 1/4"	0.6	30	thin	regular	dark-blue	100/5000
1310324	23G	1"	0.6	25	thin	regular	dark-blue	100/5000
1310325	24G	1"	0.55	25	regular	regular	purple	100/5000
1310326	25G	1"	0.5	25	regular	regular	orange	100/5000
1310327	27G	3/4"	0.4	19	regular	regular	grey	100/5000
1310328	19G	1"	1.1	25	regular	regular	cream	100/5000
1310329	19G	1 1/2"	1.1	40	thin	regular	cream	100/5000
1310330	19G	2"	1.1	50	thin	regular	cream	100/5000
1310331	18G	1 1/2"	1.2	40	thin	regular	pink	100/5000
1310332	18G	2"	1.2	50	regular	short	pink	100/5000
1310333	30G	1/2"	0.3	13	regular	regular	yellow	100/5000
1310334	27G	1/2"	0.4	13	regular	regular	grey	100/5000
1310335	26G	3/8"	0.45	10	regular	intradermal	brown	100/5000
1310336	26G	1/2"	0.45	13	regular	regular	brown	100/5000
1310337	26G	5/8"	0.45	16	thin	regular	brown	100/5000
1310338	22G	1"	0.7	25	thin	regular	black	100/5000
1310339	22G	1 1/2"	0.7	40	thin	regular	black	100/5000
1310340	22G	2"	0.7	50	regular	regular	black	100/5000
1310341	21G	5/8"	0.8	16	thin	regular	dark-green	100/5000
1310342	21G	1"	0.8	25	thin	regular	dark-green	100/5000
1310343	21G	2"	0.8	50	thin	regular	dark-green	100/5000
1310344	20G	1"	0.9	25	thin	regular	yellow	100/5000
1310345	16G	1 1/2"	1.6	40	thin	regular	white	100/5000
1310350	19G	1 1/2"	1.1	40	thin	short	cream	100/5000
1310351	18G	1 1/2"	1.2	40	thin	short	pink	100/5000
1310352	16G	1 1/2"	1.6	40	thin	short	white	100/5000
1310353	15G	1 1/2"	1.8	40	regular	short	blue gray	100/5000
1310354	14G	1 1/2"	2.1	40	regular	short	light-green	100/5000